**NOTICE D'UTILISATION APPLICATEUR SPIRE'IT®**  
**DIRECTIONS FOR USE SPIRE'IT® APPLIER**  
**INSTRUCCIONES DE USO APLICADOR SPIRE'IT®**  
**INSTRUÇÕES DE UTILIZAÇÃO DO APLICADOR SPIRE'IT®**  
**GEBRAUCHSANWEISUNG APPLIKATOR SPIRE'IT®**  
**ISTRUZIONI PER L'USO APPLICATEORE SPIRE'IT®**

# Spire'it®

Applicateur Spire'it®  
 Spire'it® applier  
 Aplicador Spire'it®  
 Aplicador Spire'it®  
 Applikator Spire'it®  
 Applicatore Spire'it®

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**TO BE CAREFULLY READ AND STRICTLY OBSERVED** by the surgeon and the staff of the operating theatre.

This instruction for use is not a reference document on surgical stapling or suturing but an information about the use of Spire'it® suture applier. The use of this device has to be carried out exclusively in operating rooms with requested aseptic conditions mentioned in this directions for use. The use of this device is strictly reserved to specialized surgeons. The non-observance of these recommendations may contribute to a reduced life duration of the device. The surgeon is responsible for any complication that may result from erroneous indication, incorrect operative technical or inadequate aseptic precautions. The patient must be informed of the precautions that he will have to take after the implantation as well as the contra-indications and undesirable side-effects that may occur. Clinical results are depending on the selected material, the surgical indication and the quality of the surgical operation.

**INDICATIONS :**

MICROVAL's Spire'it® suture applier is used in general endoscopic or open surgery for meshes fixation, during hernia repair and approximation of tissues in various surgical specialties.

**DESCRIPTION OF THE DEVICE :**

MICROVAL's Spire'it® suture applier is a sterile device designed for a multiple use during one same operation. Appliers for endoscopy are to be used with appropriate sized trocar sleeve diameter 5mm minimum.

MICROVAL's Spire'it® suture applier allows to implant memory shape nickel-Titanium alloys (NITINOL) having a double spire 4mm diameter shape.

MICROVAL's Spire'it® suture applier is delivered sterile ready to use, with 10 suture spires.

The applier implant one fixation at each pressure on the trigger (A). At expulsion of its place, the spire closes itself in spire suture shape, fixes or approximates tissues and/or meshes.

(A) Trigger, (B) extremity,(C) cartridge inside B, (D) rotation indicator at the rear of the applier.

**INSTRUCTIONS FOR USE :**

Check good compatibility of all instruments to each other. Remove MICROVAL's Spire'it® suture applier of its packaging according usual techniques for sterile products.

To avoid any damage, do not drop the applier on the sterile area. During endoscopic operation, introduce the applier through an adapted sleeve trocar (5mm diameter minimum).

Place the tip of the applier at the center of the operating site and perpendicular to tissues or mesh. In the case of approximation of tissues, use a grip forceps before any spire implementation.

Apply light pressure to the tissue site during Spire'it® suture placement. During placement of a spire internally, apply an external counter-pressure to the area being fixed whenpossible to ensure secure suture placement.

To fire the applier, simultaneously press the tip against the tissue or mesh and squeeze completely the trigger (A) and release to set the instrument for the next use.

Do not apply excessive force when squeezing the trigger the applier may jam. If a Spire'it® suture becomes jammed remove all the applier.

Check the fixed area to ensure the good shape and place of the Spire'it® . Repeat all the above steps along the length of the area to be approximated and/or fixed until the procedure is complete.

**SINGLE USE CARTRIDGE (C) :**

A rotation indicator (D) located at the rear of the applier allows to spot on the condition of the cartridge. Check the initial position of this indicator at the beginning of the procedure. Once this indicator made one 360° rotation, the 10 fixation spires have been implanted.

**PRECAUTIONS OF USE :**

An excessive pressure on the tip of the applier during firing may jam and block it.

After implantation of any Spire'it® suture with the applier, always ensure carefully the haemostasis. It is then possible to care the bleeding with electro-coagulation or manual suture.

**DISCARD AFTER USE**

All contaminated applier has to be decontaminated if it has to be return at MICROVAL.

**CONTRA-INDICATIONS :**

Do not use Spire'it® suture applier on vascular, nervous tissues nor on hard bodies (bones).

It is however possible to place Spire'it® sutures above underlying bones, vessels or viscera, but there must be a distance of not less 4.5 mm from the surface of the tissue to the underlying structure.

Do not use Spire'it® suture applier to perform haemostasis. Do not use Spire'it® suture applier on a site where it is impossible to visually check haemostasis.

Store at room temperature  
 Avoid prolonged exposure to elevated temperature.

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