



# NEW PANTHER SMART STAPLER

Smart reload identification  
Smart tissue recognition  
Smart firing control

Exclusively Distributed by:



Lap Tech Medical Sdn Bhd  
Tel: +603 - 7880 7179  
Website: [www.laptechmedical.com.my](http://www.laptechmedical.com.my)  
Email: [admin@laptechmedical.com.my](mailto:admin@laptechmedical.com.my)

[www.pantherhealthcare.com](http://www.pantherhealthcare.com)

## ELECTRONICALLY ASSISTED

Chip incorporated in the reloads for the most precise firing

Innovative 45 degrees  
Non-Stop articulation  
system

## SMART HANDLE

Universal handle ready  
to use with the most  
extended portfolio  
of reloads



One touch firing  
and opening system

Universal handle  
with better design  
and ergonomics  
for easy 360 degrees  
knobbing



# P3H

## TECHNOLOGY P3H

- Lateral diffusion of the fluids
- Stepped cartridge surface
- Better tissue compression

The P3H technology guarantees better stapling delivering better hemostasis and safety

## BENEFITS

- Enhanced anvil delivers better compression & manipulation of the tissue
- Stapling with P3H technology at the areas of greater tissue tension
- Hemostasis & safety



Purple Black



Brown



Curved Tip Reloads for the most challenging procedures



## FEATURES

Stepped cartridge



Drop-shaped staple groove



Progressive staple pushing technology



Enhanced anvil



Stronger clamping force



Four staples at the end of the staple line



Interchangeable cartridge



Coating





# NEW GENERATION ENDO STAPLING SYSTEM

Product Name	Product	Description	Color
Smart Stapler	PEAS	Powered Endoscopic Cutter Stapler with Non-Stop articulation system (Short Handle)	
	PEAM	Powered Endoscopic Cutter Stapler with Non-Stop articulation system (Standard Handle)	
	PEAL	Powered Endoscopic Cutter Stapler with Non-Stop articulation system (Extra long Handle)	
Articulating Reloading units	CADD-30EVS	Endo Articulating reloading unit 30 mm for Vascular/Medium tissue with compression range from 0.88 to 1.5mm	 Brown 12mm
	CADD-30ENTS	Endo Articulating reloading unit 30 mm for Medium/Thick tissue with compression range from 1.5 to 2.25mm	 Purple 12mm
	CADD-45EVS	Endo Articulating reloading unit 45 mm for Vascular/ Medium tissue with compression range from 0.88 to 1.5mm	 Brown 12mm
	CADD-45ENTS	Endo Articulating reloading unit 45 mm for Medium/Thick tissue with compression range from 1.5 to 2.25mm	 Purple 12mm
	CADD-45EPTS	Endo Articulating reloading unit 45 mm for Extra Thick tissue with compression range from 2.25 to 3 mm	 Black 15mm
	CADD-60EVS	Endo Articulating reloading unit 60 mm for Vascular/Medium tissue with compression range from 0.88 to 1.5mm	 Brown 12mm
	CADD-60ENTS	Endo Articulating reloading unit 60 mm for Medium/Thick tissue with compression range from 1.5 to 2.25mm	 Purple 12mm
	CADD-60EPTS	Endo Articulating reloading unit 60 mm for Extra Thick tissue with compression range from 2.25 to 3 mm	 Black 15mm
Articulated Reloads for endo linear cutter stapler	CADE-30EVS	Endo Articulating reloading unit 30 mm with Curve Tip and variable-sized staple height 2.0 mm, 2.5 mm, 3.0 mm	 Brown 12mm
	CADE-30ENTS	Endo Articulating reloading unit 30 mm with Curve Tip and variable-sized staple height 3.0 mm, 3.5 mm, 4.0 mm	 Purple 12mm
	CADE-45EVS	Endo Articulating reloading unit 45 mm with Curve Tip and variable-sized staple height 2.0 mm, 2.5 mm, 3.0 mm	 Brown 12mm
	CADE-45ENTS	Endo Articulating reloading unit 45 mm with Curve Tip and variable-sized staple height 3.0 mm, 3.5 mm, 4.0 mm	 Purple 12mm
	CADE-60EVS	Endo Articulating reloading unit 60 mm with Curve Tip and variable-sized staple height 2.0 mm, 2.5 mm, 3.0 mm	 Brown 12mm
	CADE-60ENTS	Endo Articulating reloading unit 60 mm with Curve Tip and variable-sized staple height 3.0 mm, 3.5 mm, 4.0 mm	 Purple 12mm
	CADE-60EPTS	Endo Articulating reloading unit 60 mm with Curve Tip and variable-sized staple height 4.0 mm, 4.5 mm, 5.0 mm	 Black 15mm

Exclusively Distributed by:



Lap Tech Medical Sdn Bhd  
 Tel: +603 - 7880 7179  
 Website: [www.laptechmedical.com.my](http://www.laptechmedical.com.my)  
 Email: [admin@laptechmedical.com.my](mailto:admin@laptechmedical.com.my)

# Endo Linear Cutter Stapler

WITH NON-STOP ARTICULATING LEVER

Ergonomic handle compatible with all Panther reloading units!

Innovative 45 degrees Non-Stop articulation system

Loading button on the body of the instrument

Incorporated grasping mode – grasping the tissue works by either pushing or pulling back the loop handle



# ENDO LINEAR CUTTER STAPLER

ERGONOMIC HANDLE COMPATIBLE WITH ALL PANTHER RELOADING UNITS!

Product	Description
CEAC-30	PANTHER ENDO Linear Cutter Stapler, standard handle with non-stop Articulation system and grasping mode
CEAC-45	PANTHER ENDO Linear Cutter Stapler, short handle with non-stop Articulation system and grasping mode
CEAC-60	PANTHER ENDO Linear Cutter Stapler, extra long handle with non-stop Articulation system and grasping mode

The information presented in this material is intended for product description and promotion only. Panther Healthcare has the relevant power of interpretation. Always refer to the product label and/or instruction for use before using any Panther Healthcare product. For further details and updates please refer to Panther Healthcare website: [www.pantherhealthcare.com](http://www.pantherhealthcare.com)