



FineRetriever™

Specimen Retrieval Bag

FineRetriever™

Specimen Retrieval Bag

- Minimize spillage and intraoperative contamination
- Isolate and contain infectious specimens



Indications

Safe, convenient removal of large tissue specimens such as

- Spleen
- Lung
- Ectopic Pregnancy
- Gallbladder
- Appendix
- Lymph nodes
- Ovaries
- Bowel specimens
- Other tissues structures

FineRetriever™ Specimen Retrieval Bag


Continuous Ring

The bag will remain open and not slide off the ring during specimen insertion

Radio-Opaque Rim

Clear and Strong Bag


allows for safe and convenient removal of tissue specimens such as Appendix, Gallbladder, Ovaries, Fibroid tumors etc.



FineRetriever™ Specimen Retrieval Bag

Bag Material Strength Comparison

LAGIS vs. Covidien Endo Catch™ and Ethicon EndoPouch




* Data on file from Covidien

FineRetriever™ Specimen Retrieval Bag

Weld Strength Comparison

Covidien Endo Catch™ Gold Device vs.
Ethicon EndoPouch Retriever™*



* Data on file from Covidien

FineRetriever™ Specimen Retrieval Bag

Capacity vs Trocar size LAGIS vs Covidien & Eithicon

Medium capacity	volume	Trocar to fit through
Covidien EndoCatch	220 ml	10 mm
Ethicon EndoPouch	254 ml	10 mm
LAGIS FineRetriever	320 ml	10 mm

FineRetriever™ Specimen Retrieval Bag

Capacity vs Trocar size LAGIS vs Covidien & Eithicon

Large capacity	volume	Trocar to fit through
Covidien EndoCatch	1000 ml	15 mm
Ethicon EndoPouch	-	-
LAGIS FineRetriever	1000 ml	12 mm

FineRetriever™ Specimen Retrieval Bag

- Continuous Ring
- Radio-Opaque Rim



- Strong, Tear Resistant
- Clear Pouch

Product Code	Descriptions	Bag Size	Volume
EBH-0507	FineRetriever 10mm	5 inch x 7 inch	320 ml
EB-0507	SimpleBag 10mm	5 inch x 7 inch	320 ml
EBH-0709	FineRetriever 12mm	7 inch x 9 inch	1000 ml
EB-0709	SimpleBag 12mm	7 inch x 9 inch	1000 ml

THANK YOU!