



ALPHALOOP[®]

PGLA (DMGL300) / PDO (DMML300)

- **Time saving**
- **Easy to Handle** (Patent Number 10-01170520)
- **Stable Position Control**

● Description

This product is absorbent medical suturing loop device attached aluminum tube. This product is a surgical instrument to be used for tying and suturing the blood vessels and tissues during endoscopic surgery.



● Features

- PGLA sutures are absorbed in 56~70 days.
- PDO sutures are absorbed in 180~210 days.
- The supporting aluminum tube provides effective position control.
- By using a screw connection, easy and stable operation.
- single-use medical device.

● Information

Code	Product Name	Ligature Size	Length (mm)	Sales Unit
DMGL300	ALPHALOOP - PGLA	USP 0	330	20/BOX
DMML300	ALPHALOOP - PDO	Metric[EP] 3.5	330	20/BOX