

# evolution<sup>®</sup>

**a new balance**  
the new lightweight surgical mesh



*Monofilament Polypropylene Surgical Mesh*

**dipromed**  
Dipro Medical Devices




Use scan app to read



[www.dipromed.eu](http://www.dipromed.eu)

Rete leggera, macroporosa, non riassorbibile in monofilamento di polipropilene con e senza strisce blu per facilitarne il posizionamento.  
Lightweight, macroporous, non-absorbable, monofilament polypropylene mesh with and without blue stripes to ease its positioning.



Confezionamento/Packaging:  
5 pezzi/scatola - 5 pieces/box  
\* BULEV5050 - BULEVB5050  
1 pezzo/scatola - 1 piece/box



### Bulev

Forma Shape	Codice Code	Misura Size
	BULEV0210	2x10 cm
	BULEV0611	6x11 cm
	BULEV0815	8x15 cm
	BULEV1515	15x15 cm
	BULEV1530	15x30 cm
	BULEV2535	25x35 cm
	BULEV2636	26x36 cm
	BULEV3030	30x30 cm
	BULEV3050	30x50 cm
	*BULEV5050	50x50 cm



### Bulev B - with blue stripes

Forma Shape	Codice Code	Misura Size
	BULEV1015B	10x15 cm
	BULEVB0611	6x11 cm
	BULEVB0815	8x15 cm
	BULEVB1515	15x15 cm
	BULEVB1520	15x20 cm
	BULEVB2025	20x25 cm
	BULEVB3030	30x30 cm
	*BULEVB5050	50x50 cm



#### Indicazioni/Indications:

Trattamento dei difetti della parete addominale in chirurgia open o laparoscopica.  
Treatment of the abdominal wall defects in open or laparoscopic surgery.

### P3 evolution

Forma Shape	Codice Code	Misura Size
	P3EM	4,5x10 cm
	P3EML	6x12 cm
	P3EM614	6x14 cm
	P3EF	4,5x10 cm
	P3EFL	6x12 cm
	P3EF614	6x14 cm


### P3 evolution B - with blue stripes

Forma Shape	Codice Code	Misura Size
	P3BEM	4,5x10 cm
	P3BEML	6x12 cm
	P3BEF	4,5x10 cm
	P3BEFL	6x12 cm

#### Indicazioni/Indications:

Rete presagomata per il trattamento delle ernie inguinali maschili (con foro) e femminili in chirurgia open.  
Pre-shaped mesh for treatment male inguinal hernia (with hole) and female in open surgery.

### TP bidimensional plug B - with blue stripes

Forma Shape	Codice Code	Misura Size
	TPBM	Ø 5 cm
	TPBL	Ø 7 cm

#### Indicazioni/Indications:

Tappo piano per la riparazione delle ernie inguinali e femorali, primarie e recidive.  
Flat plug for inguinal and femoral hernias, primary and recurrent.

Leggere attentamente le istruzioni per l'uso e le etichette per indicazioni, modalità d'impiego, precauzioni, controindicazioni, pericoli ed avvertenze.  
Carefully read the instruction for use and labels for indications, directions for use, precautions, contraindications, hazards and warnings.

Il sistema qualità per la progettazione, produzione e distribuzione è conforme alla norma UNI CEI EN ISO 13485:2012. I dispositivi recano marcatura CE rilasciata ai sensi dell' Allegato II della Direttiva Europea 93/42/CEE.

The quality system for design, production and distribution is in compliance with UNI CEI EN ISO 13485:2012. The devices are CE marked under Annex II of European Directive 93/42/EEC.

Fabbricante: Manufacturer:

  
Dipromed  
Dipro Medical Devices

DIPRO MEDICAL DEVICES s.r.l.  
www.dipromed.eu  
e.mail: info@dipromed.eu

Via Ciriè 22/a - 10099  
San Mauro Torinese  
(TO) - Italy

tel. +39 011 822 0298  
+39 011 669 3010  
fax. +39 011 822 3557